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```
In [ ]: %matplotlib ipynb
# The function of the previous line is to allow dynamic view of the plots

# Some base functionalities of python can be loaded through these two packages
import os
import re

# Pathlib enable to generate path valid for both Windows and Unix systems
# https://docs.python.org/3/library/pathlib.html
from pathlib import Path

# Numpy is a package containing several mathematic useful function,
# the value of pi, etc. as we will use it several time we will call it "np"
import numpy as np

# Load various scientific-related function (such as linear regression ...)
from scipy import stats

# Fabio is a package allowing us to read and write .edf images
import fabio

# Pandas can manage tables and databases, we will use it only to load tables
import pandas as pd

# lmfit is a fitting function package, we will use it to fit Pseudo-Voigt
import lmfit

# Matplotlib is a well-known plotting package
from matplotlib import pyplot as plt
from matplotlib import cm # Will be use for generating a colormap

# pyFAI is the package we will use for the calibration of the
# detector geometry and the data reduction (azimuthal integration)
import pyFAI
```

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```
In [ ]: # This is the base repository (In this case "Exercise_CeXS_Workshop", so
# If you want to move the "Script" directory to another location,
# along with all your other scripts for instance, this can become useful.
#path_base = Path('/Volumes', 'LaCie', 'TD_DRX_IMOHEM')
#path_base = Path ('/Users/rg248280/Desktop/Diffraction_synchrotron/TD Ec
path_base = Path ('/Users/raphie/Desktop/TD Ecole IMOHEM')

# Path to the directory containing the images from the detector (in .edf f
path_raw = Path(path_base, "raw_images")

# Path to the directory containing utilities for the data analysis
# (such as the experiment notebook, detector information, etc.)
path_ut = Path(path_base, "utilities")

# Path to the future integrated 1D curves, that we will compute from the
path_int = Path(path_base, "integrated")

# Path to the analysis results (graphs, dataset of stress-strain ...)
path_res = Path(path_base, "results")
```

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```
In [ ]: def get_wavelength(E):
    # Energy should be in keV !!!
    # Return the wavelength in m
    # -----
    # Plank Constant
    h = 6.62607004e-34 # m2 . kg / s
    # light celerity
    c = 299792458 # m / s
    # Joule -> eV
    eV = 6.242e+18 # eV / J

    # Conversion Energy -> wavelength
    wavelength = ((h*c)/(E*1e3 / eV)) # m

    return wavelength
```

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```
In [ ]: def integ_azm(IMG, Scan_path,
                    save_folder, Save_pyFAI,
                    N_pts, r_rng, azm_sections,
                    method_pyFAI = 'cython', integ_unit = '2th_deg'):

    # method_pyFAI = 'cython' --> parallel computing (on the pixel treatment)
    # method_pyFAI = 'csr' --> sequential computing
    # Most software work with 2theta so I keep it this way...

    # Paths...
    fileName = Scan_path.stem

    # Azimuthal integration
    if azm_sections == 1:
        # If there is only one azimuthal section, we use integrate1d()
        res = ai.integrate1d(IMG,
                             npt = N_pts,
                             unit = integ_unit,
                             mask = mask, # (1 for masked pixels, and 0 for not)
                             radial_range = r_rng,
                             method = method_pyFAI,
                             correctSolidAngle = True,
                             filename = save_path)

    elif azm_sections > 1:
        # Else, if there are several azimuthal section, we use integrate2d()
        res = ai.integrate2d(IMG,
                             npt_rad = N_pts,
                             unit = integ_unit,
                             mask = mask, # (1 for masked pixels, and 0 for not)
                             radial_range = r_rng,
                             npt_azim = azm_sections,
                             method = method_pyFAI,
                             correctSolidAngle = True,
                             filename = save_path)

    # Probably return more will be useful
    return (res.radial, res.intensity)

print("*****")
print("Integraton function defined")
print("-----")
```

```
In [ ]: # Number of integration points
# -----
N_pts = 2000

# Loading the image file of the sample
# -----
fileName = "CR1bnonirr_0"

path_IMG = Path(path_raw, fileName + ".edf")
IMG = fabio.open(path_IMG)
IMG = IMG.data

# Azimuthal integration
# -----
save_path = Path(path_int, fileName + ".dat") # It should be .dat files (

print(f"Integrating {path_IMG} ...")
res = ai.integrate1d(IMG,
                    npt = N_pts,
                    unit = "2th_deg", # We want to use 2 theta in degree
                    mask = mask, # (1 for masked pixels, and 0 for valid
                    filename = str(save_path))
print(f"saved at {save_path}")
```

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Exemple with a dummy `for:` loop:

```
for ii in range(0,10):
    print(ii)
print("The loop is over now")
```

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```
res = ai.integrate1d(IMG,
                    npt = N_pts, # Number of
integration points
                    unit = "2th_deg", # We want
to use 2 theta in degree
                    mask = mask, # previously
loaded (1 for masked pixels, and 0 for valid pixels
is the pyFAI standard)
                    azimuth_range = delta_rng, #
the range of the azimuthal angle (define at the
beginning of each loop)
                    radial_range = r_rng, # the
radial range
                    filename = save_path) # The
name of the file to save the data

print(f"saved at {save_path}")
```

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```

In [ ]: # Radial range
# -----
r_rng = [19.5, 21] # deg

# Azimuthal step
# -----
delta_step = 2 # deg
# There is no need to go from 0 to 360° as the detector do not cover
# the whole diffraction ring (measurement was in reflection).
# We will then integrate within these bounds:
delta_bounds = [210, 232]
# And then we will have N_delta steps:
N_delta = int(np.abs((delta_bounds[1] - delta_bounds[0]) / delta_step))

# Determination of the number of integration points
# -----
pix_per_pts = 1

pix_size = ai.pixell # pixel size of the detector (m)
D = ai.dist # Distance between the sample and the detector (m)
r_lenght = r_rng[1] - r_rng[0] # deg

##### --- A COMPLETER --- #####
Dtheta_pix = # size of a pixel in terms of 2theta (deg)
N_pts = # The number of integration points to perform

# Loading the image file of the sample
# -----

##### --- A COMPLETER --- #####

# Azimuthal integration
# -----

##### --- A COMPLETER (vous pouvez utiliser la partie TIP du dessus) ---

```

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```
In [ ]: # Define a range of 2theta to crop the data around the [1 1 0]

##### --- A COMPLETER --- #####

# Plot the cropped [110] peak

##### --- A COMPLETER --- #####
```

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```
In [ ]: # File to load
fileName = "CR1bnonirr_0"
# Azimuthal angles to load (first bound)
azimuth_toLoad = np.arange(delta_bounds[0], delta_bounds[1], delta_step)

# We need to initialize "I" with np.zeros to the correct shape
I = np.zeros((len(Tth), len(azimuth_toLoad)))

# Then load all the file in a for loop
# and store them into I (without forgetting to apply the crop)
for ii in range(0, len(azimuth_toLoad)):
    path_file = Path(path_int, fileName + "_" + str(azimuth_toLoad[ii]) +
    print(f"loading {path_file} ...")
    DATA = pd.read_csv(path_file, skiprows=22, delimiter = " ", skipiniti

    I[:,ii] = DATA.values[crop,1]

print("All done!")

# We can also print them all now, in a for loop again:
cmap = cm.get_cmap('jet', len(azimuth_toLoad))

fig = plt.figure()
for ii in range(0, azimuth_toLoad.shape[0]):
    plt.plot(Tth, I[:,ii], marker="o", linestyle="-", alpha=0.4, markersize=10,
             color=cmap(ii), label=f"[{azimuth_toLoad[ii]}, {azimuth_toLoad[ii]}]")
plt.yscale("log")
plt.xlabel(r"$2\theta$ (degree)")
plt.ylabel("Intensity (arbitrary unit)")
plt.legend()
fig.show()
```

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```
In [ ]: # We will call this model "peak_model"
peak_model = lmfit.models.PseudoVoigtModel(prefix='peak_')
```

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```

In [ ]: # This is not a build-in function of lmfit (at my knowledge), so we need
def background(x, y0, y1, y2):
    return(y0 + y1 * x + y2 * x**2)
# I called the three parameters of this functions "y_0", "y_1" and "y_2"

# We now define it as a lmfit model:
bkg = lmfit.Model(background)

# In order to fit at the same time the background and the peak, we need to
model = bkg + peak_model

# This time we will update the parameters "pars" with the parameters of the
# (if we had generated the parameters using .make_params(), the Pseudo-Voigt
pars.update(bkg.make_params())
pars['y0'].set(value = 62, min = 1e-08, max = 1000)
pars['y1'].set(value = 1e-3, min = 1e-08, max = 100)
pars['y2'].set(value = 1e-3, min = 1e-08, max = 10)

I_calc_init = model.eval(pars, x=Tth)

# Plot
cmap = cm.get_cmap('jet', azimuth_toLoad.shape[0])

fig = plt.figure()
ii = 0 # We will start with ii = 0 (delta = [160, 170])
plt.plot(Tth, I[:,ii], marker="o", linestyle="none", alpha=0.4, markersize=10,
         color=cmap(ii), label=f"Exp. data [{azimuth_toLoad[ii]}, {azimuth_toLoad[ii+1]}]")

plt.plot(Tth, I_calc_init, '-', color = [40/255, 190/255, 70/255], label="Initial guess")

plt.yscale("linear")
plt.xlabel(r"$2\theta$ (degree)")
plt.ylabel("Intensity (arbitrary unit)")
plt.title("Initial guess - model with a Pseudo-Voigt profile + polynomial")
plt.legend()
#fig.show()

```

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```

In [ ]: ii = 0 # We will start with ii = 0 (delta = [160, 170])

# The intensity value to be fitted need to be selected, we will store the
yy = I[:,ii]

# A common practice for XRD peak fitting is to weight the fit by the 1/I
# This way, the fit will not accord too much importance to the high value
# (We should however be careful that none of the yy values are equal to 0
wgt = np.sqrt(1/yy)

# Fitting
out = model.fit(yy, pars, x=Tth, weights=wgt) #fitting

# The fitted intensity can be found with:
I_fit = out.best_fit

# To evaluate the goodness of the fit, the value of Rwp is often used (ag
# Here is the formula to calculate it:
Rwp = (np.sqrt(np.sum(wgt*(yy-out.best_fit)**2) / np.sum(wgt*yy**2))) *

# Plot - I propose this type of plots, allowing to see the residuals (I_e
cmap = cm.get_cmap('jet', azimuth_toLoad.shape[0])

fig = plt.figure()
grid = fig.add_gridspec(4, 1, hspace=0)

ax1 = plt.subplot(grid[0:3, 0])
plt.plot(Tth, I[:,ii], marker="o", linestyle="none", alpha=0.4, markersize=
        color=cmap(ii), label=f"[{azimuth_toLoad[ii]}, {azimuth_toLoad[ii]}]
plt.plot(Tth, I_fit, '-', color = [240/255, 70/255, 40/255], linewidth=1.
plt.plot(Tth, background(Tth, out.best_values["y0"], out.best_values["y1"]

plt.yscale("linear")
plt.xlabel(r"$2\theta$ (degree)")
plt.ylabel("Intensity (arbitrary unit)")
plt.title(f"[111] peak fit (Rwp = {np.round(Rwp*100)/100} %)")
plt.legend()

ax2 = plt.subplot(grid[3, 0], sharex=ax1)
plt.plot(Tth, yy - I_fit, 'o', markersize = 1.5, alpha = 0.6,\
        color = [30/255, 30/255, 30/255], label = "residuals")
plt.legend(loc='best')
plt.xlabel(r"$2\theta$ (deg)")
plt.ylabel("Residuals (a.u.)")

#fig.show()

```

```

In [ ]: ii = 0 # We will start with ii = 0 (delta = [160, 170])

# The intensity value to be fitted need to be selected, we will store the
yy = I[:,ii]

# A common practice for XRD peak fitting is to weight the fit by the 1/I
# This way, the fit will not accord too much importance to the high value
# (We should however be careful that none of the yy values are equal to 0
wgt = np.sqrt(1/yy)

# Fitting
out = model.fit(yy, pars, x=Tth, weights=wgt) #fitting

# The fitted intensity can be found with:
I_fit = out.best_fit

# To evaluate the goodness of the fit, the value of Rwp is often used (ag
# Here is the formula to calculate it:
Rwp = (np.sqrt(np.sum(wgt*(yy-out.best_fit)**2) / np.sum(wgt*yy**2))) *

# Plot - I propose this type of plots, allowing to see the residuals (I_e
cmap = cm.get_cmap('jet', azimuth_toLoad.shape[0])

fig = plt.figure()
grid = fig.add_gridspec(4, 1, hspace=0)

ax1 = plt.subplot(grid[0:3, 0])
plt.plot(Tth, I[:,ii], marker="o", linestyle="none", alpha=0.4, markersize=
        color=cmap(ii), label=f"[{azimuth_toLoad[ii]}, {azimuth_toLoad[ii]
plt.plot(Tth, I_fit, '-', color = [240/255, 70/255, 40/255], linewidth=1.
plt.plot(Tth, background(Tth, out.best_values["y0"], out.best_values["y1"

plt.yscale("linear")
plt.xlabel(r"$2\theta$ (degree)")
plt.ylabel("Intensity (arbitrary unit)")
plt.title(f"[111] peak fit (Rwp = {np.round(Rwp*100)/100} %)")
plt.legend()

ax2 = plt.subplot(grid[3, 0], sharex=ax1)
plt.plot(Tth, yy - I_fit, 'o', markersize = 1.5, alpha = 0.6,\
        color = [30/255, 30/255, 30/255], label = "residuals")
plt.legend(loc='best')
plt.xlabel(r"$2\theta$ (deg)")
plt.ylabel("Residuals (a.u.)")

#fig.show()

```

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```
In [ ]: print(out.fit_report(show_correl=False))
```

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```
In [ ]: E = 17.6 # keV
wavelength = get_wavelength(E) * 1e10 # Ang

##### --- A COMPLETER --- #####
d_110 = # Ang # d-spacing values measured for [1 1 0] plane

# We can also calculate the center of all the delta_rng
delta_center = azimuth_toLoad + delta_step/2 # deg

# Make a plot of the d_110 spacing versus the delta_center

##### --- A COMPLETER --- #####
```

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```
In [ ]: epsilon = np.log(d_110/d0_110) # strain (without unit)

# Make the plot of the strain epsilon versus sin2(psi):
##### --- A COMPLETER --- #####
```

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```

In [ ]: # We will use the package scipy, from wich we load only "stats":
slope, intercept, r_value, _, std_err = stats.linregress(sin2_psi, epsilon)
R2 = r_value**2

# Young modulus
E = 280 # GPa

# Poisson ratio
v = 0.22

##### --- A COMPLETER --- #####
# We can then calculate the stress:
sigma =

# We can estimate an error, based on the error on the slope of the linear
# However, this error do not take into account the uncertaininties of the
# and other error sources that propagate so the errorr is minimized...
sigma_err = (std_err) * E / (1+v)

print("-----")
print(f"Sigma_11 = {round(sigma*100)/100} GPa    (+/- {round(sigma_err*100)/100} GPa")
print("-----")

fig = plt.figure()
plt.plot(sin2_psi, epsilon*100, 'o', color="r", linestyle="none", alpha=0.5)
plt.plot([0, 10], (intercept + np.multiply(slope, [0, 10]))*100, 'r-')
plt.xlim([0, 1])
plt.ylim([0, 1])
plt.xlabel(r'$\sin^2(\Psi)$')
plt.ylabel(r'$\epsilon_{[111]}$ (%)')
plt.title("Strain of the [110] plane of Cr - Regression")
plt.gca().set_aspect('equal')

```

That's all folks